

+ +

# Dehydrated Normal Human Plasma

**D** This standard Army-Navy plasma package contains the equipment to infuse 250cc of plasma within five minutes. Only trained medical officers should administer plasma. Plasma must reconstitute for five minutes before use. Rehydrated plasma must be used within three

---

Contents of This Package:

Sterile Water.....250cc      Warning: Ensure the filter is in place  
Dehydrated Plasma.....250cc      between plasma and injection needle  
Transfer Needle, Air Filter, Injection Needle      before using.

Property of the U.S. Government

**PARKE, DAVIS AND CO.**  
**DETROIT, MICHIGAN**

+ +

# Dehydrated Normal Human Plasma

**D** This standard Army-Navy plasma package contains the equipment to infuse 250cc of plasma within five minutes. Only trained medical officers should administer plasma. Plasma must reconstitute for five minutes before use. Rehydrated plasma must be used within three

---

Contents of This Package:

Sterile Water.....250cc      Warning: Ensure the filter is in place  
Dehydrated Plasma.....250cc      between plasma and injection needle  
Transfer Needle, Air Filter, Injection Needle      before using.

Property of the U.S. Government

**PARKE, DAVIS AND CO.**  
**DETROIT, MICHIGAN**